VERESS NEEDLE

AGUJA DE VERESS

AIGUILLE DE VERESS

VERESS-NADEL

AGO DI VERESS

AGULHA DE VERESS

VERESSNAALD

VERESS-NÅL

VERESS İĞNE

BEAONH VERESS

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VERESS NEEDLE

	DO NOT REUSE
STERILE EO	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
RX ONLY	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION

Please read all warnings and instructions before use.

DESCRIPTION

The Veress Needle is a sterile and single-use product. It incorporates a spring-loaded blunt stylet mechanism similar to the needle. It is used to establish peritoneum prior to trocar and cannula insertion in laparoscopic procedures. The Veress Needle, available in 120mm and 150mm lengths, has applications in gynecological laparoscopy and other laparoscopic procedures.

CONTRAINDICATIONS

This device is not intended for use when endoscopic techniques are contraindicated.

INTENDED USE

The Veress Needle consists of a stainless steel 14 gauge needle and spring-loaded stylet that retracts backward while the needle is pushed through the abdominal tissue. The stylet then advances forward once the end of the needle exits through anterior wall of the abdomen. The Veress Needle has a stopcock with a luer lock connector for insufflation tubing attachment.

WARNINGS AND PRECAUTIONS

- Federal (U.S.A) law restricts this device to sale, distribution, and use by or on the order of a physician.
- This device is EO sterilized and for single use only. Please do not re-use or re-process to avoid irreversible injury such as cross-infection or injury caused by device malfunction.
- Procedures for endoscopic surgery should be performed only by qualified and trained physicians familiar with endoscopic surgical techniques.
- A thorough understanding of techniques and principles involved in electrosurgical procedures is necessary to avoid burn or shock hazards to patient and/or operator.
- Do not use when the package is opened or damaged.
- Store this device in a dry, clean and safe place. Keep the device away from flammable and combustible chemicals or materials.
- Please follow universal guidelines upon medical waste disposal to avoid any potential cross-infection, environmental contamination or injury caused by sharps waste.

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INSTRUCTIONS FOR USE

- Prior to insertion, inspect the Veress Needle to ensure that the stylet is moving freely in and out of the needle with minimal pressure applied to the stylet. Do not attempt to use the Veress Needle if the stylet does not return over the sharp end when tension on the stylet is halted.
- The Veress Needle is inserted into the abdomen per standard insertion technique.
- Check to ensure that the Veress Needle is in the abdominal cavity and insufflate to desired pressure.
- Once insufflation is complete and the desired abdominal pressure is accomplished, remove the Veress Needle and proceed with procedure. Please follow instructions provided by the manufacturer of insufflator.
- Upon completion of the endoscopic procedure, dispose of this device in accordance with local regulations.

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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©2021 All Rights Reserved.
DeRoyal, the DeRoyal logo, and Improving Care. Improving Business. are trademarks or registered trademarks of DeRoyal Industries, Inc.