## KITNER LAPAROSCOPIC DISSECTION SPONGE

	<b>DO NOT REUSE</b>
STERILE/EO	<b>STERILIZED WITH ETHYLENE OXIDE</b>
	<b>NOT MADE WITH NATURAL RUBBER LATEX</b>
RX ONLY	<b>FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.</b>

### IMPORTANT INFORMATION

Please read all instructions and warnings before use. Correct application is essential for proper functioning of the product.

### INTENDED USE

The DeRoyal® Kittner Laparoscopic Dissection Sponge is an x-ray detectable sponge intended to be inserted through a trocar by a licensed practitioner for tissue dissection, mobilization, fluid absorption, moisture retention and tissue protection during various surgical procedures.

### ⚠️ WARNINGS

- Product should not be used if sterile packaging is damaged or opened. If damage is found, remove the product from the operative area to prevent unintended use, and return to manufacturer.
- This device is for single use only. It is not intended for disinfection and/or subsequent re-use, which can result in microbial contamination causing health deterioration of the patient.
- This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure, which could result in patient illness, injury or death.
- Do not use if this device is damaged or contaminated.
- If sponge cannot be retrieved through trocar, remove trocar and retrieve through incision site.
- Failure to remove sponge from the patient may result in a foreign body reaction.
- Count all devices before and after the procedure prior to surgical closure. In the event a device cannot be located, an x-ray can be used to locate the devices. Only the radiopaque markers are visible on imaging. The size and position of the radiopaque markers may impact their visibility.
- Avoid cutting the sponge because fragments without x-ray detectable material may enter the surgical site. Fragments left in the surgical site may result in an unintended adverse reaction.
- This device shall be used in accordance with physician instructions, hospital protocols and this IFU. Read all sections of the IFU prior to use. Improper use of these devices may cause serious injury.

### DIRECTIONS

- UPON OPENING PRODUCT PACKAGES** – Immediately verify count of sponges when package is opened.
- Insert sponge through trocar and use for tissue dissection, mobilization, fluid absorption, moisture retention and tissue protection. Remove sponges when necessary and replace.
- PRIOR TO CLOSING WOUND**– Verify sponge count is the same as it was prior to procedure.

### STORAGE AND TRANSPORT CONDITIONS

Do not store product at extreme temperatures or in a moist/damp environment; doing so may damage the product which could cause a device malfunction and/or injury to the patient.

	<b>KEEP DRY</b>
	<b>KEEP AWAY FROM SUNLIGHT</b>

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

### WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

Manufactured by:  
DeRoyal Industries, Inc.  
200 DeBusk Lane  
Powell, TN 37849 USA  
888.938.7828 or (601) 865.938.7828  
www.deroyal.com

EMERGÓ EUROPE  
Prinsessengracht 20  
2514 AP The Hague  
The Netherlands

Part#0-2020G  
Revised 2/2021

©2021 DeRoyal Industries, Inc.

All Rights Reserved. DeRoyal, the DeRoyal logo, and Improving Care. Improving Business. are trademarks or registered trademarks of DeRoyal Industries, Inc.

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