



## COVADERM® GTL

silicone composite wound dressing

DeRoyal has created the ideal secondary dressing with our COVADERM® family of wound products. Save time and money by replacing gauze, tape, and scissors with our innovative non-adherent wound dressings.

With rounded corners and more shapes and sizes than ever, COVADERM® **GTL** with a silicone adhesion layer are here to support you in creating an ideal wound healing environment.



## COVADERM® GTL

silicone composite wound dressing

- Multi-layered silicone composite wound dressing
- Dressings may be removed without causing damage to the wound or to the surrounding skin
- Ideal for use in traditional wound applications
- Rounded edges to reduce "roll-up"
- Open mesh design is coated with a soft silicone adhesive to increase absorption and reduce wound trauma
- Replace gauze, tape, and scissors with an easier solution
- Conforms easily to wound

## COVADERM® GTL


silicone composite wound dressing



† Conformable adhesive-coated silicone tape

‡ Soft pad for drainage absorption

† Protective silicone contact layer



Product #	Size	Qty/Bx	Qty/Cs	HCPC
46-401S	4" x 4"	10	50	A6203
46-402S	6" x 6"	10	50	A6203
46-403S	6" x 8"	10	50	A6204
46-411S	2.5" x 2.5"	10	50	A6203
46-412S	4" x 5"	10	50	A6203
46-408S	4" Diameter w/2" Radial Slit	25	200	A6203

### indications:

Minor abrasions, lacerations, scalds, and burns	I.V. sites, tubular sites, incisions
Dermal ulcers	More severe lacerations, cuts and abrasions
Second degree burns	Other superficial wounds
Donor sites	

### contraindications:

Third degree burns	Individuals who are sensitivity/allergy to silicone adhesive
Heavily bleeding wounds	

## application

1. Peel back edge of folded release paper.
2. Center the contact layer (gauze) over the wound.
3. Anchor exposed edge of tape onto skin, then carefully peel off rest of release paper.
4. Smooth all tape edges onto skin.



DeRoyal®, DeRoyal logo, Improving Care. Improving Business.®, and Covaderm® are registered trademarks of DeRoyal Industries, Inc.

\*All study and clinical data claims are on file.

Copyright© 2018, DeRoyal. All rights reserved.

DeRoyal Industries, Inc. | 200 DeBusk Lane, Powell, TN 37849 USA

800.DEROYAL or 865.938.7828 | www.deroyal.com

Reprint # 0-2425 | Rev. 1/18

