



STEP 1
Replace the existing OR mattress with the non-skid table mattress pad designed for the Steep-T System.



STEP 2
Place the patient lift pad and head support pillow into correct position to receive the patient. Align the perineal cutout of the lift pad with the cut out of the mattress table pad.

NOTE: The bottom 6 inches of the lift pad is not covered with the non-skid backing.



STEP 3
Position the shorter of the two straps across the lower half of the pad and attach to both sides of bed rail.

NOTE: Tighten the straps to compress the two pads together. If patient needs to be moved using the patient pad, release the hook and loop table strap before moving them.



STEP 4
Place the lift sheet over the table strap to align with the patient's lower lumbar region. The lift sheet can be used to reposition the patient and secure the patient's arms.

NOTE: The patient lift pad can be used to reposition the patient if needed. Re-secure the table strap once the patient is in final position.



STEP 5
Apply arm protectors by placing the patient's arm on the pad with the strap toward the patient's body and secure using the hook and loop closure. Repeat on the other arm. Secure both arms by wrapping the lift sheet over the patient's padded arms and tucking under the OR table pad.



STEP 6
Place the chest strap across the patient and secure to the bed railing. Final assessment of the patient position should be agreed upon by the team before prepping and draping.



STEP 7
The patient is now secure for table tilt and angulation.