
EN: HAND CONTROL MITTEN

MANUFACTURED BY:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (001) 865.938.7828
www.deryoyal.com

MedEnvoy Global B.V.
Prinses Margrietplantsoen 33 – Suite 123
2595 AM The Hague
The Netherlands

UK Responsible Person
MedEnvoy UK Limited
85, Great Portland Street, First Floor
London, W1W 7LT
United Kingdom

PART#0-716G I REVISED 10/23

©2023 DeRoyal Industries, Inc.
All Rights Reserved. DeRoyal and the DeRoyal logo
are registered trademarks of DeRoyal Industries, Inc.

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FR: MOUFLÉ DE CONTENTION

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IT: MANOPOLA DI CONTROLLO MANO

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NL: BESCHERMWANT

SE: MITSAR FÖR HANDKONTROLL

TR: EL KONTROL ELDİVENİ

EL: MONOKOMMATO PERIORIÓSTIKO GANTI

HAND CONTROL MITTEN

	SINGLE PATIENT USE
	NON-Sterile
	MEDICAL DEVICE
	NOT MADE WITH NATURAL RUBBER LATEX
RX ONLY	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION

Please read all instructions, warnings and cautions before use. Correct application is essential for proper functioning of the product. Use only on the person it was provided to by a healthcare professional and only for the use it was intended.

INTENDED USE

The DeRoyal® Hand Control Mitt is intended to help prevent patients from disturbing medical treatment.

CONTRAINDICATIONS

- Never use on a patient if an I.V. or wound site could be compromised by the device.
- Do not use on a patient who is or may become highly aggressive or agitated.

WARNINGS

- Improper application or use of this device may result in serious injury or death.
- Check to ensure that mitts are properly secured and circulation is not compromised.
- Check to ensure mitts are not damaged, missing components, or contaminated, and that all hook and loop and/or buckle closures function properly.
- Use hospital approved bed safety guidelines. If securing the device to a bed or chair, always secure straps to the portion of the bed or chair that moves with the patient. Do not secure straps to bed side rails.
- Always keep excess straps out of patient's reach. Check the patient regularly per hospital's policies and procedures.

CAUTIONS

- This product is to be fitted a properly licensed practitioner who is familiar with the purpose for which it is intended. This practitioner is responsible for providing wearing instructions and cautions to other healthcare practitioners or healthcare providers involved in the patient's care.
- The practitioner must be trained and demonstrate competency in using this device per the instructions provided and per facility policies, state and federal regulations. The practitioner should perform a

complete assessment of the patient to ensure the device is necessary. Use a restraint only when all other alternative methods of treatment and intervention have failed to meet the patient's needs.

- Be sure to follow your facility's policies and guidelines for frequency of patient monitoring and documentation.
- Consult your prescriber immediately if you experience sensation changes, unusual reactions, swelling or increased pain while using this brace. Discontinue use if pressure injuries develop.

INSTRUCTIONS FOR USE

- Release the hook strap and open the wrist closure.
- Place the patient's hand in the mitt with the padded side touching the patient's palm. If there are finger separators, make sure fingers are placed into the appropriate slots.

3. Secure wrist closure, with either hook and loop strap or laces. If using hook and loop strap, pull the strap until the wrist closure fits snug around the patient's wrist. Secure the hook to the loop side of the strap. If using laces, gently pull both laces until the wrist closure fits snug around the patient's wrist, and tie the laces as if tying a shoelace. **NOTE:** Do not over-tighten, two fingers should fit between the device and the patient to ensure circulation is not compromised.

- IF APPLICABLE:** The optional bed attachment strap(s) may be used for further security. Feed the strap(s) through the loop hole on the wrist closure. Feed the end of the strap through the sea wave buckle and pull completely through until the teeth engage the strap. Secure the strap(s) using quick release ties to the portion of the bed frame that moves with the patient.

QUICK RELEASE TIES

Wrap the strap(s) around the anchor once. Take the loose end of the strap(s) and create a loop. Cross the loop over the other strap(s) that is wrapped around the anchor and feed it through opening as if tying a shoelace. Hold that loop and create another loop with the loose strap(s) and feed it through the first loop. Pull the final loop and the strap(s) attached to the cuff in opposite directions to form a quick-release tie.

CLEANING AND/OR MAINTENANCE (IF APPLICABLE)

Please refer to hospital approved guidelines for cleaning and disinfecting of the device.

STORAGE AND TRANSPORT CONDITIONS

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROVAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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