## PRO-GLIDE™ KNEE ORTHOSIS

ORTESIS DE RODILLA PRO-GLIDE™

ORTHESE DE GENOU PRO-GLIDE™

KNIEORTHESE PRO-GLIDE™

ORTOSI PER GINOCCHIO PRO-GLIDE™

ÓRTESE PARA JOELHO PRO-GLIDE™

PRO-GLIDE™ KNIEORTHESE

PRO-GLIDE™ KNĀORTOS

PRO-GLIDE™ DİZ ORTEZİ

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Manufactured by:  
DeRoyal Industries, Inc.  
200 DeBak Lane  
Powell, TN 37849 USA  
615.938.7020 or (801) 865.938.7228  
www.derooyal.com

EMERGO EUROPE  
Prinsesgracht 20  
2514 AP The Hague  
The Netherlands

Part #1-4105G

Revised 2/21

©2021 DeRoyal Industries, Inc.

All Rights Reserved. DeRoyal, the DeRoyal logo, and Improving Care. Improving Business. are trademarks or registered trademarks of DeRoyal Industries, Inc.

## PRO-GLIDE™ KNEE ORTHOSIS

	<b>SINGLE PATIENT USE</b>
	<b>NON-STERILE</b>
	<b>NOT MADE WITH NATURAL RUBBER LATEX</b>
<b>RX ONLY</b>	<b>FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.</b>

### IMPORTANT INFORMATION

Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury or re-injury inherent with the use of any brace. Use only on the person it was provided by a healthcare professional and only for the use it was intended.

### INTENDED USE

This knee orthosis was designed for the treatment of joint stiffness and established soft tissue contractures which may have resulted from trauma or immobilization.

### WARNING

- Inspect brace for damaged or missing components before use.
- Discontinue use and consult your prescriber if the brace or its components break or become damaged.
- Prolonged usage can result in muscle atrophy in the user. The user may require (1) a weaning period from the brace to build muscle strength (2) regular physical therapy exercises to maintain the strength of the leg (3) a limitation on duration of use, or (4) some combination of these. The prescriber should communicate these risks to the user and determine the treatment plan.

### CAUTIONS

- A physician or properly licensed practitioner (a "prescriber") who is familiar with the use and purpose of this brace must fit it to the user. The prescriber has a duty to provide wearing instructions and risks related to the use of this brace to other healthcare practitioners treating the users and the users themselves, including duration of use. The instructions provided in this sheet do not supersede hospital protocol or direct orders of the prescriber. Use only as directed.
- Inspect the fit of the brace on the patient to ensure fit is proper. **DO NOT OVERTIGHTEN.** Check regularly to ensure circulation is not compromised.
- Consult your prescriber immediately if you experience sensation changes, unusual reactions, swelling or increased pain while using this brace. Discontinue use if pressure injuries develop.
- While not an absolute contraindication to the use of this brace, take special care if the user is diabetic or has poor circulation. These users may have decreased skin sensitivity and are at greater risk for poor peripheral circulation and pressure injuries.
- Once the sizing and range of motion limitations are set, use caution while wearing to ensure these settings are not inadvertently adjusted.
- Metal components could cause injury to the user or others. Exercise care in situations where contact is likely, such as contact sports.

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### INSTRUCTIONS FOR USE

- To adjust proximal cuff, lay the orthosis open with all anterior straps unfastened and the posterior straps loosened. Turn adjustment screw clockwise to lower the distal struts below the point of contracture. Apply orthosis to the posterior side of the leg and align hinge axis with the axis of the knee. Firmly fasten proximal cuff closure above the knee. Secure corresponding hook and loop straps over proximal cuff closure.
- To adjust anterior thigh pad, tuck pad between thigh and proximal cuff and secure it with the strap. Next tighten the posterior thigh strap. **NOTE:** Optimal fit requires the body of the hinge to rest securely along the longitudinal axis of the femur. The center bracket can be adjusted upward or downward as required using 3/8" nut driver.
- To adjust the distal cuff for extension, firmly fasten around the calf. Align the distal struts slightly posterior to the longitudinal axis of the tibia. Fasten both posterior calf straps. Apply the center strap over the anterior side of the leg.
- To adjust proximal cuff for flexion, the distal cuffs and straps must be removed and reversed with the hoop and loop closure to the anterior and straps over the posterior.
- To adjust the distal cuff for flexion, firmly fasten distal cuff around calf. Align distal struts slightly anterior to longitudinal axis of the tibia. Adjust both anterior calf straps. Apply the center strap over the posterior side of the leg.
- To adjust tension for extension and flexion, increase or decrease force by turning the adjustment screw. Turn the adjustment screw to achieve the desired force level. Make sure that the tension setting is the same on both hinges.
- To convert the tension from extension to flexion, turn the force adjustment screw and reverse the tension on the cable. Continue turning screw until no slack is present in the cable and the desired tension has been reached. **NOTE:** When the cable is wrapped around the bottom portion of the pulley, the orthosis is in flexion. When wrapped around the top portion of the pulley, the orthosis is in extension.
- The locking feature enables the patient to don and doff the orthosis without reducing the tension setting. Remove the lock pin from its storage holder and place in the locking hole (3 available). Gradually move arm against the applied tension force while pushing inward on the lock pin. The pin will slide into a receiving hole with a click when it is aligned making a "click" sound. Repeat for other side. To remove the lock pin, move the distal struts against the applied tension force to take pressure off the pin. By simultaneously pulling and rotating the lock pin, it can be removed and placed back into its holder.

### THESE INSTRUCTIONS APPLY TO PRODUCT NUMBERS: 4105 SERIES

### CARE INSTRUCTIONS

Clean hardware with a nonabrasive cleaner. Straps and pads may be hand washed in mild detergent in lukewarm water. Air dry.

### PARTS

Air Foam Cushion Anterior Thigh Pad Proximal Cuff Center Strap Distal Cuff Center Bracket Adjustment Tool Tension Adjustment Tool Tension

### Window

### FORCE ADJUSTMENTS

Force measured at 7" from axis. Force readings are approximate.

### FORCE INDEX FORCE

	lbs	metric
1	6 lbs	2.7 kg
2	10 lbs	4.5kg
3	13.5 lbs	6.12 kg
4	17 lbs	7.7 kg
5	20 lbs	9 kg
6	54 lbs	10.88 kg

### STORAGE AND TRANSPORT CONDITIONS

	<b>KEEP DRY</b>
	<b>KEEP AWAY FROM SUNLIGHT</b>

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

### WARRANTY

DeRoyal products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

### Window

### FORCE ADJUSTMENTS

Force measured at 7" from axis. Force readings are approximate.

### FORCE INDEX FORCE

	lbs	metric
1	6 lbs	2.7 kg
2	10 lbs	4.5kg
3	13.5 lbs	6.12 kg
4	17 lbs	7.7 kg
5	20 lbs	9 kg
6	54 lbs	10.88 kg

### STORAGE AND TRANSPORT CONDITIONS

	<b>KEEP DRY</b>
	<b>KEEP AWAY FROM SUNLIGHT</b>

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

### WARRANTY

DeRoyal products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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Ciddi olaylar, hastanın ikamet ettiği ülkedeki yetkili makama ek olarak DeRoyal Industries, Inc.'ye bildirilmelidir.