DEROM® ELBOW DYNAMIC ORTHOSIS
ORTESIS DE CODO DINÁMICA DEROM®
ORTHÈSE DYNAMIQUE DU COUDE DEROM®
DEROM® DYNAMIC ELLBOGENORTHESE
ORTOSI DINAMICA PER GOMITO DEROM®
ÓRTESE DINÁMICA PARA COTOVELO DEROM®
DEROM® DYNAMISCHE ELLEBOOGORTHESE
DEROM® DYNAMISK ARMBÅGSSORTOS
DEROM® DİRSEK DİNAMİK ORTEZİ
ΔΥΝΑΜΙΚΗ ΟΡΘΩΣΗ ΑΤΚΩΝΑ DEROM®

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DEROM® ELBOW DYNAMIC ORTHOSIS

	SINGLE PATIENT USE
	NON-STERILE
	NOT MADE WITH NATURAL RUBBER LATEX
RX ONLY	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION
Please read all instructions, warnings and cautions before use. Correct application is essential for proper product function and to reduce the risk of injury or re-injury inherent with the use of any brace. Use only on the person it was provided to by a healthcare professional and only for the use it was intended.

INTENDED USE
The DeRoyal® Range Of Motion (DeROM®) Dynamic Elbow Orthosis was designed for the treatment of joint stiffness, and established soft tissue contractures, which may have resulted from trauma, immobilization, or neurological disorders. Protocols for use vary depending on the diagnosis and should be determined after consulting with the patient's physician.

CONTRAINDICATIONS
Not to be used if any of the following conditions exist: Severe spasticity, severe osteoporosis, phlebitis or thrombophlebitis. Diabetes requires frequent skin checks.

- CAUTIONS**
- A physician or properly licensed practitioner (a "prescriber") who is familiar with the use and purpose of this brace must fit it to the user. The prescriber has a duty to provide wearing instructions and risks related to the use of this brace to other healthcare practitioners treating the users and the users themselves, including duration of use. The instructions provided in this sheet do not supersede hospital protocol or direct orders of the prescriber.
 - Use only as directed by medical professional.
 - Prolonged usage can result in muscle atrophy in the user. The user may require (1) a weaning period from the brace to build muscle strength (2) regular physical therapy exercises to maintain the strength of the elbow (3) a limitation on duration of use or (4) some combination of these.
 - Inspect brace for damaged or missing components before use.
 - Discontinue use and consult your prescriber if the brace or its components break or become damaged.
 - Inspect the fit of the brace on the patient to ensure fit is proper. **DO NOT OVERTIGHTEN.** Check regularly to ensure circulation is not compromised.
 - Consult your prescriber immediately if you experience sensation changes, unusual reactions, swelling or increased pain while using this brace. Discontinue use if pressure injuries develop.
 - Take special care if the user is diabetic or has poor circulation as these users may have decreased skin sensitivity and are at

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

greater risk for poor peripheral circulation and pressure injuries.

INSTRUCTIONS FOR USE
IMPORTANT: The tension indicator dial (10) must be set to "0". To set the tension Indicator dial (10) push the tension release lever (6) in the direction of the arrow. At the same time turn the tension indicator dial to "0". All slack must be taken out of the cable before applying tension. After applying orthosis, push in on the black hinge knob (2) and rotate the pulley until the cable is taut.


- TO SET UP FOR EXTENSION ASSIST:** The tension indicator dial (10) must be set to "0". See **"IMPORTANT"** above. Position the red cable (1) around the bottom of the pulley. Push in and hold the black hinge knob (2). At the same time push the forearm strut (3) to its lower most position. Release the black knob, and pull the forearm cuff to its upper most position. Repeat until the red cable is wrapped around the bottom of the pulley, and tension is felt on the cable.
- TO SET UP FOR FLEXION ASSIST:** The tension indicator dial (10) must be set to "0". See **"IMPORTANT"** above. Position the red cable (1) around the top of the pulley. Push in and hold the black hinge knob (2). At the same time pull the forearm strut (3) to its upper most position. Release the black knob, and push the forearm cuff to its lower most position. Repeat until the red cable is wrapped around the top of the pulley, and tension is felt on the cable.
- DONNING THE UNIT:** Position the arm in its end range. Unfasten forearm straps (4) and biceps straps (5). Hold the arm in the neutral (thumb pointed up) position. Place orthosis on the arm with the mechanism on the lateral side (outside) of the arm, positioning the axis of the hinge centered at the elbow joint.

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
4. ADJUST THE FOREARM CUFF LENGTH: Remove the adjustment tool (6) from the rear of the plastic housing. Loosen the forearm cuff screws (7), and position the forearm cuff (9) around the forearm, and cinch the two straps (4) for a snug, comfortable fit. **NOTE:** Set tension indicator dial to "0" prior to taking slack out of the red cable. See **"IMPORTANT"**. All slack must be taken out of the cable before applying tension. Push in on the black hinge knob (2) and rotate the pulley until the cable is taut.

- MOLD THE BICEP CUFF:** Mold the bicep cuff (8) around the bicep, and cinch the two straps (5) for a snug, comfortable fit.
- MOLD THE FOREARM CUFF:** Mold the forearm cuff (9) around the forearm, and cinch the two straps (4) for a snug, comfortable fit. **NOTE:** Set tension indicator dial to "0" prior to taking slack out of the red cable. See **"IMPORTANT"**. All slack must be taken out of the cable before applying tension. Push in on the black hinge knob (2) and rotate the pulley until the cable is taut.
- APPLYING TENSION:** Remove the adjustment tool (6) and insert it into one of the holes on the edge of the tension indicator dial (10). Rotate the dial in the direction of the arrow. Each "click" increases the tension. Adjust tension to patient tolerance. Release the adjustment tool by snapping it into the rear of the plastic housing.
- RELEASING TENSION:** Push the tension release lever (6), at the back end of the plastic housing, in the direction the arrow indicates.
- TENSION LOCKING CONTROL (TLC):** Once the appropriate tension level is determined, unscrew the TLC knob (11), and position it in the threaded hole nearest the tension indicator arrow, to prevent applying excessive tension. Screw the TLC knob completely in.
- SETTING ROM LIMITERS:** Pull the ROM stops (12) outward, and out of the stop hole. Reposition to the desired position, and insert into the appropriate stop hole.
- REMOVING ORTHOSIS:** Release tension per step #8. Open the two forearm straps (4), and the two biceps straps (5) by pulling on the colored end of each strap.

CARE INSTRUCTIONS
Clean foam cuffs and exterior of hinge with non-abrasive cleaner. Straps may be machine washed. Air dry prior to wearing.

STORAGE AND TRANSPORT CONDITIONS	
	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY
DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROWAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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All Rights Reserved. DeRoyal, the DeRoyal logo,
and DeRom are trademarks or registered
trademarks of DeRoyal Industries, Inc.