

NEGATIVE PRESSURE WOUND THERAPY DRESSING KITS

MD	MEDICAL DEVICE
	DO NOT REUSE
	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION
Please read all warnings and instructions before use.
DeRoyal intends these products to be used for patients that would benefit from Negative Pressure Wound Therapy (NPWT).





INDICATIONS
Negative Pressure Wound Therapy Dressing Kits are to be used in conjunction with an NPWT pump to promote wound healing by removal of wound exudate, debris, and infectious material. NPWT Dressing Kits may be used during surgery or at the patient's bedside and are indicated for home use under a clinician's care.

CONTRAINDICATIONS

- Necrotic tissue with eschar present
- Unexplored or non-enteric fistulas
- Untreated osteomyelitis
- Wounds containing malignant tissue
- Exposed arteries, nerves (including vagus nerve), blood vessels, veins, or internal organs
- Exposed anastomotic site

WARNING

- Read and understand all NPWT pump and kit warnings, cautions and directions completely and carefully before use. Failure to do so may result in serious injury to either the patient or the user.
- These kits are intended to be used by licensed practitioners familiar with NPWT. User should have a thorough understanding of wound care protocols and procedures to determine appropriate application and treatment and removal.
- Do not reuse, reprocess, or re-sterilize any part of the kit. This kit is for single use only. Reuse, reprocessing, or re-sterilization may damage the kit, may lead to pump failure or contamination, and may result in user or patient injury, patient infection, or patient death.
- Do not use if package is opened or package or kit is damaged. Handle with care to avoid damage to the kit packaging and contents during transportation and storage.
- Thoroughly inspect the pump, including its accessories and cables, for breaks, cracks or other damage or defects prior to use.
- If any part of the package or kit is damaged, replace the entire kit.
- These kits are designed for use with DeRoyal® Negative Pressure Therapy Units. Use with units from other manufacturers could lead to serious injury.
- Place dressing to minimize friction.
- Ensure exudate flow is not obstructed (tubing is kinked or stenosis in the tubing).
- Do not place Transreal® over the top of the dome. This may cause a closed system.
- Secure tubing of suction dome to minimize fall hazard.

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CAUTIONS

- Precautions should be taken for patients who are or have:
 - Difficult wound hemostasis
 - On anticoagulation or active bleeding
 - Close proximity of blood vessels, organs, muscle, wounds in close proximity and fascia requiring adequate protection. These areas should be completely covered and protected to avoid direct contact with wound dressing.
 - Irradiated infected vessels and tissue
 - Bony fragments
 - Untreated malnutrition
 - Non-compliant or non-responsive
 - Paralysis associated with the affected area or spinal cord injuries
 - Circumferential dressing application
 - Use near vagus nerve
- Infected wounds should be closely monitored and may require more frequent dressing changes than non-infected wounds. View exudates for any sign of infection, worsening infection, or other complications.
- If dressing adherence occurs, take proper steps to avoid damaging granulating tissue and minimizing pain.
- If infants, children, some small adults and elderly patients should be closely monitored for excessive fluid loss and dehydration. In addition, patients with heavily exuding wounds or large wounds should be monitored carefully because they may have a risk of excessive fluid loss and dehydration.
- Protect fragile periwound skin with adhesive sealing films, hydrocolloid, or transparent film dressing. Do not cover intact skin with gauze/foam dressing to avoid risk of maceration or injury to tissue.
- Patients requiring:
 - APH
 - Defibrillation
 - Hyperbaric chamber

DIRECTIONS FOR USE DURING THE PROCEDURE

1. Read the **WARNINGS** and **CAUTIONS** for all devices thoroughly and complete the steps described prior to the procedure.
2. Inspect the packaging for damage or defects. Replace the entire kit if the packaging is damaged.
3. Turn pump off prior to initiating dressing change.
4. Follow your facility's protocol regarding wound site preparation and debridement. Cleanse wound gently and carefully with saline until the wound is free of debris and dead tissue. Gently blot dry the wound edges, periwound margin, and all skin to be covered by transparent dressing. **WARNING:** Do not use if package is opened or package or kit is damaged. Handle with care to avoid damage to the kit packaging and contents during transportation and storage.
5. Follow your facility's protocol regarding wound site preparation and debridement. Cleanse wound gently and carefully with saline until the wound is free of debris and dead tissue. Gently blot dry the wound edges, periwound margin, and all skin to be covered by transparent dressing. **WARNING:** Do not use if package is opened or package or kit is damaged. Handle with care to avoid damage to the kit packaging and contents during transportation and storage.
6. Size and cut the DeRoyal® Transreal® transparent dressing to cover the primary dressing, making sure the primary dressing is completely covered and the Transreal® dressing overlaps onto the intact skin by at least one inch on all edges of the wound. Proper Transreal® dressing application is the key to experiencing fewer alarms and seal leakage issues.

7. Gently press the top of the Transreal® dressing to ensure adhesive adheres to the intact skin. Remove the top carrier layer of the Transreal® dressing.
8. Gently pinch and lift the Transreal® dressing over the center of the wound filler. Cut a 2 to 2.5 cm hole in the Transreal® dressing using either scissors or a scalpel.
9. Do not cut a slit as a slit could collapse and block the suction when pressure is applied.
10. Remove the first release liner backing from the Topdraw® dome. Center the Topdraw® dome over the hole in the Transreal® dressing and press the adhesive of the dome liner down onto the Transreal® dressing.
11. Remove the second release liner from the Topdraw® dome and press gently to ensure the Topdraw® dome adhesive is attached securely to the Transreal® dressing.
12. Remove the clear carrier film from the top of the Topdraw® dome.
13. Pat around the Topdraw® dome to recheck the seal is completely adhered to the Transreal® dressing.
14. Connect the Topdraw® dome silicone tubing to suction canister connector.
15. Turn on the NPWT system, making sure that the pressure is set at the proper setting per the clinician's order.
16. Once the NPWT system is turned on, the primary or secondary dressing, black foam, white foam, and/or gauze, Topdraw® dome pad, and Transreal® dressing will start to compress. A compressed look means the Transreal® dressing is properly adhered to the skin, the Topdraw® dome is properly adhered to the Transreal® dressing, and negative pressure is reaching the wound. **WARNING:** If dressing is not compressed, check for leaks in the system. **WARNING:** Check the wound dressing at shift change and every 8 hours to ensure negative pressure is reaching the wound.
17. **WARNING:** Document the type and number of wound dressing pieces that are placed in the wound.
18. Change dressing according to physician's orders.
19. After use, discard the entire kit according to your facility's protocol.

STORAGE AND TRANSPORT CONDITIONS
Store at room temperature and avoid excessive heat and humidity.

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY
DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **NO WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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