SURGICLICK® LIGHT HANDLE COVER

	DO NOT REUSE
	MEDICAL DEVICE
	DO NOT USE IF PACKAGE IS DAMAGED
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.
STERILE PRODUCTS ONLY	
	STERILIZED WITH ETHYLENE OXIDE
	DO NOT RESTERILIZE
	SINGLE STERILE BARRIER SYSTEM

IMPORTANT INFORMATION

Please read all instructions, warnings, and cautions before use. Correct application is essential for proper functioning of product.

INTENDED USE

The product is intended to cover a surgical light handle and maintain a sterile field during medical or surgical procedures.

INTENDED USERS

The product is intended to be used only by a licensed practitioner or trained healthcare professionals.

PATIENT TARGET GROUP(S)

Device is non-patient contacting.

INDICATIONS

Use this product whenever surgical lights are required during a surgical case and will need to be adjusted by personnel who are operating within the sterile field.

WARNINGS

- Ensure cover fits light handle securely. A loose cover could fall during a procedure and compromise the sterile field.
- DO NOT** reuse, reprocess, or re-sterilize. This device is for single use only. Reuse, reprocessing, or re-sterilization may damage the device, and may result in user or patient injury, patient infection, or patient death.
- Product should not be used if sterile packaging is damaged or opened. If damage is found, remove the product from the operative area to prevent unintended use.
- DO NOT** use damaged product or product with missing components.

INSTRUCTION FOR USE

- Use only with SurgiMate® Adapting Handle.
- Open package using aseptic technique.
- Slide two fingers into the handle portion of the light handle cover to create the correct shape and to enable you to slide the cover over the light handle.
- Place sterile cover over the shaft of the SurgiMate® Adapting Handle.
- With a quick, upward motion, snap cover over metal ring on Adapting Handle. Listen for "click" during application.
- At the end of the procedure, remove the light handle cover from the light handle by pulling down on the flange portion of the light handle cover.
- Dispose of properly in accordance with facility protocol after use.

STERIZATION INFORMATION

NON-STERILE PRODUCT: DeRoyal intends that non-sterile product for use in sterile environments be further processed by our customers, including further packaging and/or sterilization according to their own validated processes. These products are to be sterilized by ethylene oxide only, as other sterilization methods have not been validated by DeRoyal and may damage the product which could cause a device malfunction and/or injury to the patient.

DISPOSAL

Dispose of any used device according to local, state, and federal laws and regulations. For safe disposal of devices, follow your facility's protocol.

STORAGE AND TRANSPORT CONDITIONS

Do not store at extreme temperatures or in moist/damp environment; doing so may damage the product which could cause a device malfunction and/or injury to the patient.

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal products are warranted for one hundred and twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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