- FOR USE WITH YSI 400 SERIES ELECTRONICS ONLY
- PARA USO SOLO CON COMPONENTES ELECTRÓNICOS SERIE YSI 400
- À N'UTILISER QU'AVEC LE MATERIEL ÉLECTRIQUE SÉRIE YSI 400
- AUSSCHLIESSLICH ZUR VERWENDUNG MIT ELEKTRONIKGERÄTEN DER REIHE YSI 400
- DA UTILIZZARE ESCLUSIVAMENTE CON L'ELETTRONICA DELLA SERIE YSI 400
- PARA USO SOMENTE COM APARELHOS ELETRÔNICOS SÉRIE YSI 400
- UITSLUITEND VOOR GEBRUIK IN COMBINATIE MET ELEKTRONICA VAN SERIE YSI 400
- ENDAST FÖR ANVÄNDNING MED YSI 400-SERIENS ELEKTRONIK
- SADECE YSI 400 SERİSİ ELEKTRONİK CIHAZLARLA KULLANILABİLİR
- ΑΠΟΚΛΕΙΣΤΙΚΑ ΓΙΑ ΧΡΗΣΗ ΜΕ ΗΛΕΚΤΡΟΝΙΚΕΣ ΣΥΣΚΕΥΕΣ ΣΕΙΡΑΣ YSI 400

NEONATAL TEMPERATURE PROBE

	DO NOT REUSE
	MEDICAL DEVICE
	NON-STERILE
	NOT MADE WITH NATURAL RUBBER LATEX
	RX ONLY FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION

Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury.

INTENDED USE
The DeRoyal® Neonatal Skin Surface Temperature Probe is a non-sterile device intended for connection to an infant radiant warmer to monitor the neonatal patient's skin surface temperature.

WARNINGS
Inaccurate temperature readings may result if:

- Device is damaged
- Directions for Use are not followed
- Incorrect placement or application of device
- Use of product outside of the intended use is prohibited.
- Failure of adhesive pad during use may result in incorrect temperature readings. Do not continue use until the adhesive pad is replaced.
- Discontinue use and replace device immediately if temperature readings cease during use.
- While adverse reactions are rare, the following have been associated with use of the skin temperature sensor during use or immediately after use: electrical burns and skin irritation or inflammation.
- The equipment cannot differentiate between an increase in core temperature with a cold skin (fever) and a low core and skin temperature (hypothermia). The patient's temperature should be monitored.

CAUTIONS
Inspect device prior to use for damage.
Discontinue use if an allergic reaction is observed from use of the adhesive pad cover.
Consider a less aggressive adhesive pad cover if skin irritation results upon removal of the cover.
Electric shock hazard. Do not use if wire insulation is damaged or bare wires are exposed.
Take care to ensure that the cable and connector do not get wet and proper techniques are used during electro surgical procedures to reduce radio frequency

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interference current and potential burn to the patient. An adequate electrosurgical dispersive ground electrode close to the active surgical site should be properly connected.

- DO NOT** intertwine the cables, especially the monitor cables with the electrosurgical unit's cables.
- The operation of the infant radiant warmer may be temporarily affected during electrosurgical activations. Unusual temperature readings should be checked.
- The adhesive cover on this probe should never be applied to a skin surface which has been damaged by abrasion, laceration, incision, or burning.
- Do not modify the device as it may result in damage and/or injury to the patient.
- Environmental conditions (e.g., air movement) can affect the thermal balance of the infant.

MR COMPATIBILITY

The Neonatal Skin Temperature Sensor has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Neonatal Skin Temperature Sensor in the MR environment is unknown. Scanning a patient who has this medical device may result in patient injury.

DIRECTIONS FOR USE

- Prepare the skin for adhesive pad by removing body oils and hair.
- Remove the Skin Temperature Sensor from its package.
- Peel off release liner and apply sensor with probe cover to patient's skin.
- Follow clinician's discretion for probe placement. Best practice is on the patient's stomach or side.
- Follow the directions for use of the infant radiant warmer.
- When temperature monitoring is completed, carefully remove adhesive cover from patient's skin.
- Disconnect the probe (at the connector) and discard the probe. To disconnect, grasp both connectors firmly and pull. Do not pull on cable or wire.

CLEANING INSTRUCTIONS FOR INTERFACE CABLES

- If necessary, the cable may be cleaned by wiping it using a Sani-Cloth® Plus cleaning cloth, Sani-Cloth® Bleach cleaning cloth, Sani-Cloth® AF3 cleaning cloth, or an equivalent cloth utilizing the one of the following active ingredients: n-alkyl dimethyl ethylbenzyl ammonium chlorides, n-alkyl dimethyl benzyl ammonium chlorides, and sodium hypochlorite.
- DO NOT** immerse the cable in liquids for cleaning.
- DO NOT** autoclave the cable or gas sterilize if temperature will exceed 60°C (140°F).

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OPERATING SPECIFICATIONS

- Rated output range: 25 to 45° Celsius
- The reference body site is the skin surface temperature.
- The measuring site is the skin surface.
- Accurate to ±0.2° Celsius

STORAGE AND TRANSPORT CONDITIONS

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT
	DO NOT USE IF PACKAGE IS DAMAGED
	RELATIVE HUMIDITY: 85% NON-CONDENSING
	OPERATING TEMPERATURE: -25°C TO +45°C
	STORAGE TEMPERATURE: -25°C TO +55°C

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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