ALGIDEX AG® I.V. PATCH

PARCHE I.V. ALGIDEX AG®

ALGIDEX AG® I.V. PATCH

PATCH I.V. ALGIDEX AG®

Manufactured by:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (001) 865.938.7828
www.deroyal.com

PART# 660906
Revised 6/2021

©2021 DeRoyal Industries, Inc.
All Rights Reserved. DeRoyal, the DeRoyal logo,
Improving Care. Improving Business., and Algidex
Ag are trademarks or registered trademarks of
DeRoyal Industries, Inc.

ALGIDEX AG® I.V. PATCH

	USE BY DATE
	DO NOT REUSE
	STERILIZED USING IRRADIATION
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION
Please read all instructions, warnings and cautions before use. Correct application is essential for proper functioning of the product. Use only on the person it was provided to by a healthcare professional and only for the use it was intended.

DEVICE DESCRIPTION
Algidex Ag® I.V. Patch is a sterile patch of polyurethane foam coated with an ionic silver alginate and maltodextrin matrix. Algidex Ag® I.V. Patch contains silver that is evenly distributed in the ionic state. The dressing provides a sustained antimicrobial activity in the dressing for up to seven (7) days against a broad spectrum of microbial organisms.

INTENDED USE
Algidex Ag® I.V. Patch is intended to manage exudating wounds, including abrasions and lacerations, dermal wounds, donor sites, first- and second-degree burns, surgical incisions, vascular access sites, pressure ulcers stage I-IV, stasis ulcers, and venous ulcers.

CONTRAINDICATIONS

- Algidex Ag® I.V. Patch should **NOT** be used:
 - on patients with known sensitivity to alginates or on patients with sensitivity to silver or silver compounds

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- on third degree burns
- ulcers resulting from infections
- lesions in patients with active vasculitis
- DO NOT** surgically implant.

WARNINGS

- Local wound care cannot overcome the deficits of unrelieved malnutrition, pressure or trauma, or compromised blood flow. Interventions must be instituted to reduce/relieve these factors that contribute to impaired wound healing in order to achieve the optimum outcome. When visible signs of infection are present, use these dressings in conjunction with proper systemic medical treatment to address the cause of infection.
- DO NOT** use when contraindicated.
- Discontinue use if signs of sensitivity or irritation appear.
- Product is sterile unless packaging is opened or damaged. **DO NOT** use if product is expired or if packaging is open or damaged.
- Discard excess product after initial opening.
- DO NOT** use for longer than seven (7) days.
- This product is prescription only. **DO NOT** provide it to the patient without providing complete use instructions.
- Remove product and cleanse wound thoroughly prior to hyperbaric oxygenation or any other therapy where the presence of metals is contraindicated.

CAUTIONS

- DO NOT** store in extreme temperatures or in a damp environment. This may render the product ineffective.

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INSTRUCTIONS FOR USE
ALGIDEX AG® I.V. PATCH FOAM

- May be applied immediately following initial catheter placement or during routine catheter dressing change.
- Apply Algidex Ag® I.V. Patch with silver matrix touching the catheter and skin.
- Secure patch in place with retention dressing such as gauze, transparent film, or tape.
- Patch may be worn up to seven (7) days.

REMOVAL OF ALGIDEX AG® I.V. PATCH

- Remove retention dressing, if applicable.
- Lift Algidex Ag® I.V. Patch from the catheter and discard according to institutional policy.

EXPIRATION INFORMATION

- Store Algidex Ag® I.V. Patch in a cool, dry place.
- Contents sterile unless package has been opened or damaged.
- See package label for expiration date.

STORAGE AND TRANSPORT CONDITIONS

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT
	DO NOT USE IF PACKAGE IS DAMAGED

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY
DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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INSTRUCTIONS D'UTILISATION PATCH I.V. EN MOUSSE ALGIDEX AG®
RETRAIT DU PATCH I.V. ALGIDEX AG®