DEROYAL® MANIFOLDS AND STOPCOCKS

	DO NOT REUSE
	MEDICAL DEVICE
	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

INTENDED USE
DeRoyal® manifolds and stopcocks are intended to be used during cardiac catheterization procedures for intra-arterial and intravenous administration of contrast, saline or radiographic contrast media. The product is intended for use by a properly licensed practitioner.

PRODUCT	MAXIMUM PRESSURE RATING* (PSI)	MAXIMUM PRESSURE RATING (KPA)
Low Pressure Manifolds	250 PSI	1723 kPa
Medium Pressure Manifolds	550 PSI	3792 kPa
High Pressure Manifolds (designed for procedures requiring small catheters and power injection)	800 PSI	5515 kPa

! WARNINGS

- All connections should be hand tightened, without over-tightening, to ensure their security and prevent damage to the device.
- Check all connections to ensure they are tight and secure to prevent air from entering the system.
- Examine the products carefully for entrapped air and fully de-bubble prior to injection.
- Contents supplied non-sterile should be sterilized prior to use. Sterilization should be by Ethylene Oxide according to the manufacturer's validated sterilization parameters.
- Connections should be made only with devices compliant with ISO 80369-7 for luer connections.

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- Low Pressure Manifolds are not intended for power injection as this may result in leaking and/or allow air bubbles into the system.
- The PSI ratings listed are for injections made through the end port. The side port(s) are not rated and should not be used for high pressure injections.
- Inspect packaging for any damage and product for any damage, contaminate or missing parts prior to use.

PRECAUTIONS

- DO NOT reuse. This device is not intended for disinfection and subsequent re-use, which may result in device failure or create the risk of contamination.
- The device has not been evaluated for re-processing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure which could result in patient illness, injury or death.

STERILIZATION INFORMATION

DeRoyal intends that non-sterile products for use in sterile environments be further processed by our customers, including further packaging and/or sterilization. These products are to be sterilized by ethylene oxide, utilizing the following parameters: Temperature: 120°F - 135°F (48.9°C - 57.2°C), Ethylene Oxide Gas Concentration: 625 - 800 mg/L, Exposure Time: 4 hrs minimum (Preconditioning: Temperature 100°F - 120°F (37.8°C - 48.9°C), Humidity 55-75%, Time 18hrs minimum), (Aeration: Temperature 100°F - 135°F (37.8°C - 42.2°C), Time 18hrs minimum). Other sterilization methods have not been validated by DeRoyal and may damage the product which could cause a device malfunction and/or injury to the patient.

*PLEASE NOTE: All pressure ratings are based on DeRoyal's sterilization cycle. Other sterilization may affect product performance.

STORAGE AND TRANSPORT CONDITIONS

Do not store product at extreme temperatures or in a moist/damp environment.

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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