



OA double upright KNEE BRACE

COMFORT for the way you **LIVE**

CONTROL for the way you **PLAY**

CONFIDENCE in how you **RECOVER**

 **DeRoyal**[®]

Improving Care. Improving Business.™

OA double upright KNEE BRACE

FEATURES & BENEFITS

- Features a unique thigh cuff that is flexible yet durable; the three strapping system enables optimal suspension, protection and comfort by moving with the muscles throughout range of motion.
- Thigh cuff accommodates changes due to muscle atrophy, hypertrophy, weight gain or loss, and can be adjusted to accommodate larger patients.
- Strong and lightweight fiber-reinforced composite frame.
- Interchangeable condyle pads for OA unloading, increased comfort, stability and pain control.
- Numbered strapping system for intuitive donning and doffing.
- Silicone gripped calf strap liner to help prevent migration.
- Optional flexion/extension stops for range of motion control



INDICATIONS: Mild to moderate unicompartmental osteoarthritis

SIZE	MEASUREMENT 6" ABOVE MID-PATELLA	MEASUREMENT MID-PATELLA	MEASUREMENT 6" BELOW MID-PATELLA	LEFT	RIGHT
SMALL	15½" - 18½"	11" - 13"	12¼" - 13¾"	14515105	14515205
MEDIUM	18½" - 21"	13" - 15"	13¼" - 15"	14515106	14515206
LARGE	21" - 23½"	15" - 17"	14¼" - 15¾"	14515107	14515207
X-LARGE	23½" - 25"	17" - 19"	15" - 17"	14515108	14515208
XX-LARGE	25" - 28"	19" - 22"	17" - 19"	14515109	14515209



Improving Care. Improving Business.®

DeRoyal Industries, Inc.
 200 DeBusk Ln | Powell, TN 37849
 888.938.7828 | www.deroyal.com
 Revision 6/14 | Reprint # 0-2344

