



FUNCTIONAL ACL

KNEE BRACE

LOW CONTACT | CONTACT

COMFORT for the way you **LIVE**

CONTROL for the way you **PLAY**

CONFIDENCE in how you **RECOVER**

 DeRoyal®

Improving Care. Improving Business.™

FUNCTIONAL ACL KNEE BRACE

FEATURES & BENEFITS

- Features a unique thigh cuff that is flexible yet durable and enables optimal suspension, protection and comfort by moving with the muscles throughout the range of motion (ideal for LOW-CONTACT sports/use only.)
- Thigh cuff accommodates changes due to muscle atrophy, hypertrophy, weight gain or loss, and can be adjusted to accommodate larger thighs.
- Aerospace grade aluminum frame is strong yet customizable.
- Also comes in full-frame design (ideal for CONTACT sports/use.)
- Posterior calf cuff reduces tibial pressure.
- Interchangeable condyle pads for increased comfort, stability, and protection.
- Numbered strapping system for intuitive application and removal.
- Silicone gripped calf strap liner to help prevent migration.
- Optional flexion/extension stops for range of motion control.



full frame

INDICATIONS: Mild to moderate ACL, PCL, MCL, and LCL sprains | Combined ligament instabilities

SIZE	MEASUREMENT 6" ABOVE MID-PATELLA	MEASUREMENT MID-PATELLA	MEASUREMENT 6" BELOW MID-PATELLA	LOW CONTACT LEFT	LOW CONTACT RIGHT	CONTACT LEFT	CONTACT RIGHT
SMALL	14½" - 17"	11½" - 13"	11½" - 13"	14533105	14533205	14534105	14534205
MEDIUM	17" - 19½"	13" - 15"	13" - 14½"	14533106	14533206	14534106	14534206
LARGE	19½" - 22½"	15" - 17"	14½" - 16½"	14533107	14533207	14534107	14534207
X-LARGE	22½" - 25½"	17" - 19"	16½" - 18½"	14533108	14533208	14534108	14534208
XX-LARGE	25½" - 30"	19" - 22"	18½" - 22"	14533109	14533209	14534109	14534209



Improving Care. Improving Business.®

DeRoyal Industries, Inc.
200 DeBusk Ln | Powell, TN 37849
888.938.7828 | www.deroyal.com
Revision 6/14 | Reprint # 0-2345

