




- Easy to Use
- Prevents Patient Movement
- Enhances Airway Management
- Reduces Head/Neck Angle

A unique, intuitive positioning solution to enhance airway management for patients in surgical procedures. Prevents patients from migrating/slipping on the surgical table while in Trendelenburg position up to 40° angulation.



Intuitive Design

- Reduces OR set up time (less components)
- Less risk of incorrect setup that can lead to patient injury

Deep Head Cutout

- Stabilizes the patient's head
- Maintains ideal position for airway maintenance

Positioning Kit

- Packaged with one patient chest strap
- Optional lift sheet and arm protector components available

Ramp

- Reduces head/neck angle
- Sniffing position to assist in intubation of difficult airways

Ramp Wings

- Provides extra support for shoulders

Securement Bar

- Prevents patient movement when in steep angulation

Ordering Information

Product No.	Description	Quantity
M60-AP500	Disposable Pillow and Chest Strap	2/Cs
M60-AP501	Disposable Pillow, Chest Strap, Arm Protectors, and Lift Sheet	2/Cs
M60-AP510	Securement Bars	2/Cs