


VAULT EX – How to assign a product



Badge in



Select the prescribed product




1. Close the Vault door and press “DONE” on the touch screen
2. Take tablet and proceed to patient’s room
3. Fit the product to the patient



Login to the tablet using your credentials



1. Select your patient
2. Select the physician who prescribed the product
3. Select the correct product to assign
4. Select “Dispense”



1. Select the “Signer” and then provide the tablet for their signature
2. Select “Done signing”
3. Verify green check mark appears before logging out
4. Return tablet to its designated area

Questions? Call **1-800-993-9012**.
If after hours, leave a message and your call will be returned