

GAUZE ROLLS & SPONGES

**FLUFTEX™**

FLUFTEX™



- Pre-washed fluff dried gauze
- Ideal for management of all draining wounds
- Ideal for all prepping, cleansing, and bandaging needs

 DeRoyal®



## GAUZE ROLLS & SPONGES

# FLUFTEX™

### FEATURES & BENEFITS

- Pre-washed, fluff dried gauze with superior absorption and softness for packing, prepping, cleansing, and bandaging
- Crinkled to give ideal loft for the management of all draining wounds
- Extra bulk for added protection to sensitive wound sites
- Minimal linting, with finished or folded edges to prevent thread separation
- Wide variety of sizes
- Available with Raypaque™ element
- Use as primary or secondary dressing retention



Fluftex gauze rolls and sponges are available in a wide assortment of sizes and packaging.



Fluftex's finished edges help prevent pieces of raw gauze from getting into the wound site.

<b>Fluftex Sponges</b> Stretched, diagonal measure			
Product #	Size	Qty/Bx	Qty/Cs
11-001	6" x 6.75", Medium Sterile 2's in soft pouch	40	480
11-002	7.75" x 8.75", Large Sterile 2's in soft pouch	40	400
11-003	6" x 6.75", Medium Sterile 5's in soft pouch	50	600
11-004	6" x 6.75", Medium Sterile 10's in plastic tray	10	480
11-005	6" x 6.75", Medium Non-sterile bulk	100	600
11-006	7.75" x 8.75", Large Non-sterile bulk	100	400
<b>Fluftex Rolls</b> Stretched, 6-Ply			
Product #	Size	Qty/Cs	
11-007	4.5" x 4.1 yds, Sterile in soft pouch	100	
11-008	3.4" x 3.6 yds, Sterile in soft pouch	96	
11-009	2.25" x 3 yds, Sterile in soft pouch	96	
11-010	4.5" x 4.1 yds, Sterile in plastic tray	60	
11-011	4.5" x 4.1 yds, Non-sterile bulk	100	
11-012	3.4" x 3.6 yds, Non-sterile bulk	96	
<b>Sterile with Raypaque™</b> Stretched, 6-Ply			
Product #	Size	Qty/Cs	
11-020	4.5" x 4.1 yds, Sterile in soft pouch	100	
11-021	3.4" x 3.6 yds, Sterile in soft pouch	96	
11-022	2.25" x 3.0 yds, Sterile in soft pouch	96	

TO PLACE AN ORDER OR FOR MORE INFORMATION

PLEASE CALL: **800-251-9864**

888.938.7828 (WITHIN USA)

865.938.7828 (OUTSIDE USA)

800.543.2182 (FAX)



200 DeBusk Lane  
Powell, Tennessee 37849 U.S.A.  
[www.deroyal.com](http://www.deroyal.com)

Reprint #0-969 • ©2005 DeRoyal.  
All Rights Reserved.  
Rev. 5/2005