



STRETCH NET™ TUBULAR ELASTIC DRESSING

GENERAL Stretch Net is designed to hold dressing securely in place without the use of tape. Ideal for patients who are allergic to tape or who require frequent dressing changes, Stretch Net is packaged in a variety of widths and two lengths for many applications. Avoid tight, constrictive application. Choose the correct size.

LATEX SAFE STRETCH NET

REGULAR STRETCH NET

REF. 10 YARDS STRETCHED	REF. 25 YARDS STRETCHED	REF. 10 YARDS STRETCHED	REF. 25 YARDS STRETCHED	SIZE	APPROXIMATE CIRCUMFERENCE STRETCHED	SUGGESTED APPLICATION
10-7401	10-7421	10-7101	10-7125	1	3–8"	Fingers, Toes
10-7402	10-7422	10-7102	10-7126	2	8–11"	Fingers, Wrist, Small Hand
10-7403	10-7423	10-7103	10-7127	3	11–16"	Hand, Elbow, Foot
10-7404	10-7424	10-7104	10-7128	4	16–22"	Hand, Elbow, Foot, Knee
10-7405	10-7425	10-7105	10-7129	5	22–26"	Calf, Thigh, Head
10-7406	10-7426	10-7106	10-7130	5½	26–35"	Thigh, Head, Shoulder
10-7407	10-7427	10-7107	10-7131	6	35–42"	Small Chest, Axilla, Abdomen, Thigh, Shoulder
10-7408	10-7428	10-7108	10-7132	7	42–50"	Chest, Axilla, Abdomen
10-7409	10-7429	10-7109	10-7133	8	50–58"	Chest, Axilla, Abdomen
10-7410	10-7430	10-7110	10-7134	9	58–64"	Chest, Axilla, Abdomen
10-7411	10-7431	10-7111	10-7135	10	64–72"	Chest, Axilla, Abdomen
10-7412	10-7432	10-7112	10-7136	11	72–80"	Extra Large Torso
10-7413	10-7433	10-7113	10-7137	12	80–86"	Extra Large Torso
10-7414	10-7434	10-7114	10-7138	13	86–92"	Extra Large Torso



APÓSITO ELÁSTICO TUBULAR STRETCH NET™

CONSIDERACIONES GENERALES Stretch Net está diseñado para mantener el apósito fijo en su lugar de aplicación sin recurrir al uso de cinta adhesiva. Es ideal para pacientes alérgicos a la cinta o que requieren un cambio frecuente de apósitos. Stretch Net se ofrece en diferentes anchuras y dos longitudes para varios tipos de aplicaciones. Evite aplicarlo de forma ceñida. Asegúrese de seleccionar el tamaño correcto.

STRETCH NET SIN LÁTEX

STRETCH NET COMÚN

REF. 9 M ESTIRADO	REF. 23 M ESTIRADO	REF. 9 M ESTIRADO	REF. 23 M ESTIRADO	TAMAÑO	APROXIMADA CIRCUNFERENCIA ESTIRADO	SUGERIDA APLICACIÓN
10-7401	10-7421	10-7101	10-7125	1	7,6–20 cm	Dedos de la mano y de los pies
10-7402	10-7422	10-7102	10-7126	2	20–28 cm	Dedos de la mano, muñeca, manos pequeñas
10-7403	10-7423	10-7103	10-7127	3	28–40,6 cm	Mano, codo, pie
10-7404	10-7424	10-7104	10-7128	4	40,6–55,8 cm	Mano, codo, pie, rodilla
10-7405	10-7425	10-7105	10-7129	5	55,8–66 cm	Pantorrilla, muslo, cabeza
10-7406	10-7426	10-7106	10-7130	5½	66–88,9 cm	Muslo, cabeza, hombro
10-7407	10-7427	10-7107	10-7131	6	88,9–106,7 cm	Pecho pequeño, axila, abdomen, muslo, hombro
10-7408	10-7428	10-7108	10-7132	7	106,7–127 cm	Pecho, axila, abdomen
10-7409	10-7429	10-7109	10-7133	8	127–147 cm	Pecho, axila, abdomen
10-7410	10-7430	10-7110	10-7134	9	147–162,6 cm	Pecho, axila, abdomen
10-7411	10-7431	10-7111	10-7135	10	162,6–182,8 cm	Pecho, axila, abdomen
10-7412	10-7432	10-7112	10-7136	11	182,8–203 cm	Torso extra grande
10-7413	10-7433	10-7113	10-7137	12	203–218,4 cm	Torso extra grande
10-7414	10-7434	10-7114	10-7138	13	218,4–233,6 cm	Torso extra grande



Manufactured by:

DeRoyal Industries, Inc.

200 DeBusk Lane

Powell, TN 37849 USA

888.938.7828 or (001) 865.938.7828

www.deroyal.com

PART#4-6999 REVISED 6/2023

©2023 DEROYAL INDUSTRIES, INC.
All Rights Reserved. DeRoyal, the DeRoyal logo,
and StretchNet are trademarks or registered
trademarks of DeRoyal Industries, Inc.