




## Not Pictured:

- Additional cuff panel hook and loop engagement
- Anti-migration silicone strapping provided

## Universal Recovery Post-Op Knee Brace

- Lightweight
- Enhanced locking ROM hinge makes adjustment easy, increases patient compliance, and reduces the chance of inadvertent change in ROM settings
- Easy latch-and-release lever buckles make donning and doffing easy
- Adjustable inner panel location to avoid incision sites and provide a customized fit
- Compressed foam padding provides additional cushion and increased airflow
- Universally sized with intuitive telescoping length adjustments (18" – 27")
- Numerous drop lock configurations with five convenient preset drop lock settings
- Additional cuff panel hook and loop engagement helps prevent the struts from migrating around the leg during patient movement
- Flexible panels increase patient comfort
- Soft, trimmable straps with removable y-tab for customized fit
- Extended padding protects ankle from soft tissue irritation
- Anti-migration silicone strapping provided
- XL available for 33" - 39" thigh



**Product  
Information**



**Full Application  
Video**



HCPSC CODE:  
L1832 or L1833

| Size          | Length Adjustment | Product # |
|---------------|-------------------|-----------|
| Universal     | 18" - 27"         | KB3000-01 |
| Universal, XL | 18" - 27"         | KB3000-XL |

### HCPSC DISCLAIMER

Information related to L-Codes or ICD-10 Codes for use in determining the proper Medicare or Medicaid reimbursement status for any DeRoyal product, represents the reasonable efforts of DeRoyal Industries, Inc. in determining relevant L-Codes and corresponding diagnoses for certain of DeRoyal's products. This information is provided as a no-cost service to customers of DeRoyal. The information is compiled from the descriptors and functions associated with the L-Codes and a review of ICD-10 Codes from publicly available information and anecdotal information obtained from sources familiar with the billing process as indicated above. However, only the appropriate Medicare authority can ultimately determine the appropriateness of reimbursement of a particular product, for a particular patient, given the appropriate diagnosis. DeRoyal strongly cautions the users of this information to verify the applicability of a given L-Code and/or ICD-10 Code in a particular situation prior to seeking reimbursement for a DeRoyal® product