

TRANSITION KNEE BRACE



Indications

- Post-operative stabilization and rehabilitation
- Post total knee replacement
- Post ORIF
- Fractures of the femur, tibia, fibula or patella
- Ligament reconstructive procedures
- Meniscal repairs
- Post-fracture management

Features & Benefits

- Drop-lock design offers easy access for pin adjustment
- Range of motion can be set in 10° increments
- Superior and inferior score zones for length adjustment
- Circumferential loop lock straps allow for even compression
- Trimmable softgoods available in full or cool foam options

Product Information

Product #	Description	Length Adjustment
KB5000-00	Full Foam	23" - 17"
KB5000-01	Cool Foam	23" - 17"
KB5000-02	Full Foam, X-Long	26" - 20"



200 DeBusk Lane
Powell, TN 37949
www.deroyal.com
888.938.7828 • 865.938.7828
rev 12/07